How to Make Your AED AED

READY FOR RESCUE



When seconds could mean the difference between life and death for a sudden cardiac arrest (SCA) victim, be sure your automated external defibrillator (AED) is visible, easily accessible and maintained, just like a fire extinguisher.

PLACE IT



Placement of the AED should be centralized in your facility or in an area where the most people gather. The goal is to deliver a shock to an SCA victim within 3 minutes or less from collapse. Consider installing next to or near a fire extinguisher. Never lock an AED in an office or closet.

HANGIT

REACH IT

Be sure your AED is installed at an unobstructed maximum height of 48 inches from the floor in compliance with the American Disabilities Act.

SPOTIT

Use 3-dimensional universal AED Add stickers to outer doors to show an

REGISTER IT

- With AED manufacturer; it may extend the warranty
- · With local or state EMS

CHECK IT



- Designate a person in your facility to perform regular AED status checks as indicated by the manufacturer
- Check battery and pads expiration dates for timely updating
- · Subscribe to an AED Program Management service
- Enroll in free AED Registry at National AED Registry.com

STOCK IT

Include an AED Rescue Kit with all the essentials

- ✓ Nitrile Gloves
- ✓ Barrier Mask
- ✓ Medical Scissors
- ✓ Razor
- ✓ Antiseptic Wipes
- ✓ Absorbent Towel
- 🖊 Biohazard Bag

Add AED location(s) to site maps that show fire extinguishers and exit routes, then place site maps around your facility. AEDs are designed for the public to use, so be public about it!

PROMOTE IT



USE IT

SCA leads to death in minutes if the person does not get help right away. Survival depends on people nearby calling 911, starting CPR, and using an AED as soon as possible. Having a written and well-practiced cardiac emergency

















ParentHeartWatch.org